RWD

RWD Life Science Co.,Ltd

R420 Series Veterinary Electric Turbo Ventilator/ Ventilador de turbina veterinario de la serie R420

User Manual/Manual de usuario

A

Intellectual Property Right

The intellectual property rights of this product and its user manual belong to RWD Life Science Co., Ltd (hereinafter referred to as the RWD), including, but not limited to patent, trademark, copyright etc.

RWD reserves final interpretation right of this user manual.

RWD has the right to use the instruction as confidential information. Any individual and/or organization shall not disclose the instruction of all or part of the information by any means without RWD's written permission. Nor shall any other person or organization be allowed to obtain all or part of the information of this user manual by any means.

No individual and/or organization shall publish, modify, reproduce, issue, rental, adapted, and translated into other languages without RWD's written permission.

is RWD's registered trademark or mark, these trademarks and the related security mark belong to RWD's intangible property. The use of non RWD's trademark or mark in this user manual are only for editing purposes, without other purposes, the rights belong to their respective rights owners.

Statement

RWD reserves the right to modify the content of this manual without prior notice.

RWD reserves the right to change the technology without prior notice.

RWD reserves the right to modify the product specifications without prior notice.

RWD does not guarantee the information in any forms, including (but not limited to) the responsibility of proposing the implied merchantability and suitability for a specific purpose.

RWD in only the following conditions is considered to be responsible for the safety, reliability and performance of the instruments, i.e.:

Assembly operation, expansion, adjustment, improvements and repairs carried out by RWD authorized personnel;

Relevant electrical equipment in line with national standards;

The instrument is operated according to the user manual .

RWD is not responsible for the products' safety, reliability or operation status in the following conditions:

Components are disassembled, stretched or debugged;

Non RWD authorized personnel repairs or alterations to the instruments;

Product may not in accordance with the manual.

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CONTENTS/DIRECTORIO

USER MANUAL (EN)	1
FOREWORD	1
1 - PRODUCT OVERVIEW	1
1.1 PRODUCT FEATURES.....	1
1.2 ENVIRONMENT REQUIREMENTS	2
1.3 PRODUCT PARAMETERS.....	2
1.4 INSTRUMENT COMPONENTS AND INTERFACES.....	3
1.5 PRODUCT LIST	5
1.6 INTENDED PURPOSE	5
2 - SYSTEM SECURITY	6
2.1 SAFETY INSTRUCTIONS	6
2.2 WARNINGS AND PRECAUTIONS.....	6
2.3 NOTICE FOR USE	6
2.4 SAFETY INSTRUCTIONS	7
3 - UNPACKING AND ASSEMBLY	8
3.1 UNPACKING.....	8
3.2 INSTALLATION OF THE BELLOWS COMPONENTS	8
3.3 PIPELINE CONNECTION.....	9
3.4 CONNECTION OF ACCESSORIES AND ATTACHMENTS.....	12
3.5 INSTALLATION METHOD OF EtCO ₂ WATERTRAP.....	13
4 - OPERATING INSTRUCTIONS	14
4.1 POWER-ON SELF-CHECK.....	14
4.2 MAIN INTERFACE	15
4.3 VENTILATION MODE.....	18
4.4 PARAMETER RANGE	18
4.5 USE OF THE SHUTTLE KNOB	19
4.6 PARAMETER LIST	19
4.6.1 Add/edit/delete/call parameters	20
4.7 START/END VENTILATION	21
4.8 PRESSURE-TIME DIAGRAM & CO ₂ -TIME DIAGRAM	21
4.9 ALARM SETTING	22
4.10 SETTINGS	23
4.10.1 Brightness	23
4.10.2 Language.....	23
4.10.3 Device Zero	23
4.10.4 CO ₂ Zero	24
4.10.5 CO ₂ on/off.....	24
4.10.6 Remote Pair	24
4.11 USAGE AND REMOTE CONTROL OF INSP.HOLD BUTTON	25
4.12 DATE & TIME	26

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User Manual (EN)

Foreword

First of all, thank you for choosing the R420 Series Veterinary Electric Turbo Ventilator made by RWD!

Please be sure to carefully read all the information attached before installing and using this product for the first time, which will help you use the product better.

RWD is committed to continuously improving product functions and service quality, and reserves the right to make changes to any product described in this user manual and the content of this user manual without prior notice.

For the latest product information, please contact us by telephone or mail, or visit our website (<http://www.rwdstco.com/>). Please contact RWD if you find any inconsistency between the actual product and the information contained herein during the use of the instrument, or when you have any questions or suggestions.

1 - Product overview

The R420 series veterinary electric turbo ventilator is mainly used for regular pet operation, laparoscopic surgery, thoracotomy, etc. In conjunction with anesthesia machines (can be used with other brands of anesthesia machines), it can provide animals with anesthetic gases and oxygen, and simulate and cooperate with the normal breathing of animals, so as to ensure the animals' respiratory function during the operation, and avoid the operation accident caused by the animals' weakened respiratory function during the operation.

This user manual applies to the following R420 series veterinary electric turbo ventilator products made by RWD:

- R420 Veterinary Electric Turbo Ventilator
- R420-Plus Veterinary Electric Turbo Ventilator
- R420-Pro Veterinary Electric Turbo Ventilator

The product should be operated and managed by trained professionals!

This product is only for use on animals, not for humans!

1.1 Product features

- 7-inch high-definition touch screen, visualized data, convenient operation, real-time reading of monitoring data
- Three modes available: VCV mode, PCV mode, APNEA mode
- Equipped with two bellows of different sizes (namely 300ml and 1,500ml), for animals in the weight range of 0.4~150kg
- Intelligent parameter recommendation: After setting the weight, the recommended respiratory parameter value can be obtained intelligently under the flexible setting.

- Audible alarm and alarm text message prompt, which can improve the human-machine interaction between the device and the user, and reduce human errors during the use of the device.
- Provide internal rechargeable battery power
- Remote wireless control module and sidestream end-tidal CO₂ monitoring function module are optional

1.2 Environment requirements

Please prepare the instrument operating environment under the conditions listed in the table below to ensure the operation and safety of the instrument.

	Description
Operating conditions	Temperature: 5°C ~ 40°C
	Humidity: 15%~95% (non-condensing)
	Working air pressure: 70kPa~110kPa
Storage conditions	Temperature: -20°C~55°C
	Humidity: ≤95%
	Air pressure for storage and transportation: 70kPa~110kPa
Working power supply	100~240V AC, 50/60Hz, maximum rated power: 100VA
Waterproof and dustproof grade	IPX1

1.3 Product parameters

Parameter item	Description
Dimension	290×256×167mm (±10mm, excluding bellows components)
Weight (including bellows components)	8.5kg
Touch screen	7-inch screen, screen pixel: 1024×600
Battery working time	In the VCV mode, when the tidal volume is 80ml, the respiratory rate is 12, and the I:E ratio is 1:2, it can work continuously for at least 2 hours for an 8kg animal

Ventilation parameters	Setting range
Weight range	0.4~150kg
PEAK (PCV only)	5~50 cmH ₂ O
RR	2~60 bpm
I:E	1:1.0~1:4.0
VT (VCV only)	4~1500ml
P-Trig	-10 ~ -1.0 cmH ₂ O, OFF
PEEP	3~20 cmH ₂ O

1.4 Instrument components and interfaces

Figure 1-1

No.	Part	Description
1	Bellows components	Inject fresh oxygen and anesthetic gas into animal lungs (with black folding pouch inside) Bellows components: 1,500ml, 300ml
2	Display	Display the operation information and setting parameters of the device
3	External power indicator	The power indicator will light up when the device is connected to an external power supply
4	Power switch	Press and hold for 1s to turn on the device. The device will turn on and enter the self-check interface after the device prompt tone is heart. Press and hold for 3s to turn off the device.
5	Ventilator EtCO ₂ watertrap	Optional. Used for dehydration of exhaled air
6	INSP.HOLD button	Inspiration Hold button, long press to fill the chest cavity to the desired pressure or tidal volume and hold to not allow exhalation until the button is released. For details, refer to Section 4.11
7	Shuttle knob	Rotate it to switch to select a parameter, click it to confirm the parameter, and adjust the value. For the detailed operation of the shuttle knob, refer to the description in Section 4.5


Figure 1-2

No.	Part	Description
8	30mm port	Anesthetic gas outlet port, being used to exhaust the anesthetic gas in the bellows
9	22mm port	Anesthetic gas inlet port, being connected to the anesthesia machine
10	Antenna	Receive remote control signal
11	USB-B port	Software update
12	USB port	Data import and export. Only FAT32 format USB flash drive supported
13	Pressure sensor port	Connected to 22mm port
14	Battery position	The rechargeable battery can be replaced, please refer to Section 7.3 for the replacement method
15	Exhaust gas outlet	Connected to the gas filter canister
16	Driving gas output port	Connected to 17mm port
17	Fuse box	The fuse box can be pried open to replace with a new fuse. Please refer to Section 7.4 for the replacement method
18	Power cord socket	Power input
19	17mm port	Connected to the driving gas output port

1.5 Product list

Note: Due to differences in different versions of the user manual, this product list is for reference only. Please compare the delivered parts with the packing list attached with the product when receiving. If you find any inconsistencies, please contact the after-sales service personnel of RWD immediately.

Configuration	Name		Qty	Description
Standard	Mainframe (control unit)		1	Main part of the ventilator, being used to control the operation of the ventilator
Standard	Power line		1	Power supply for the mainframe
Standard	300ml bellows		1	The 300ml bellows
Standard	300-1,500ml bellows		1	The 300-1,500ml bellows
Standard	Airway pressure connection tube - Luer - 22mm connector		1	Include clear connector. Connect to the pressure sensor port of the control unit (PRESSURE TRANSDUCER)
Standard	Y-shaped connecting pipe	Bellows connector (30 to 22mm)	1	The connector is installed on the 30mm port of the bellows cover assembly, and the connector is connected to corrugated pipe
		Assisted trachea - I.D.15.0×L200mm	1	Assisted rubber tube
		T-shaped adapter	1	Adapter joint for corrugated pipe and assisted trachea
		Reducing coupling - 15mm(M)-15mm(F)/22mm(M)	1	22 to 15mm connector, connected to the assisted trachea
		Corrugated pipe - ID 22mm, L0.3m	2	Connected to the bellows connector (30 to 22mm), the driving gas output port
Standard	Corrugated pipe - ID 22mm, L1.2m		1	Connected to the EXHAUST port and the exhaust air filtering tank
Standard	Threaded pipe - ID 22mm, L1.2m		1	Connected to the anesthesia machine and ventilator
Optional	Sidestream end-tidal CO ₂ monitoring function		1	Configurable only for R420-Plus/R420-Pro models
Optional	INSP.HOLD remote control function		1	Configurable only for R420-Pro model

1.6 Intended purpose

This product can be used in pet hospitals, agricultural and forestry colleges, scientific research units and other related occasions for animal research, treatment, and surgery.

2 - System security

2.1 Safety instructions

Caution!

- The animal and respiratory equipment must be under constant observation during respiration.
- The breathing circuit connected to the anesthesia machine should be unobstructed and easy to observe, and shall not be overlaid or otherwise impair its function (eg, improper positioning of the animal).
- Please take measures to protect rubber and rubber parts from UV light and prolonged direct sunlight to avoid embrittlement of these parts.

2.2 Warnings and precautions

- The responsibility for device fault caused by improper cleaning, maintenance and operation shall be borne by the user.
- Any disassembly and assembly of the host without authorization of RWD will enable the company not to provide any quality assurance and technical maintenance service commitment for this device. For any technical problems, please contact the authorized persons or RWD for support.

Electric shock hazard!!!

It is forbidden for non-professionals to remove the cover or panel of the device and disassemble the built-in rechargeable battery.

2.3 Notice for use

Please read this user manual carefully before using this device!

This device is designed to perform the functions mentioned in this manual. Please use and operate the device according to the corresponding instructions in the manual. This device must be checked and maintained regularly to ensure its reliable operation.

If any component is partially or completely damaged, overworn, contaminated or its service life ends up for any other reason, it shall be discontinued and replaced with a component produced by the Company or an enterprise authorized by the Company.

Warning!

Do not press the quick oxygenation button on the anesthesia machine during the inspiratory phase of the respiratory cycle, which may cause severe lung rupture of animal. During the inspiratory phase, the discharge valve in the control unit is closed, so the oxygen flush flow will be added to the air supply flow of the ventilator, and the only way for the air flow at this time is the lung of experimental animal.

The device that can not work properly shall be discontinued immediately, and the user shall not make any kind of alteration to the device. The user shall be solely responsible for any injury to the subject animals or any damage to the property (including the device itself) caused by any

failure to operate in accordance with the maintenance instructions authorized by the Company, or caused by unauthorized repair or improvement of the device and its accessories, or caused by the use of damaged accessories or unauthorized accessories.

2.4 Safety instructions

This device is only for use by the personnel with knowledge of respiratory anesthesia or the personnel trained by the manufacturer.

The device shall be installed in a place away from water source, well ventilated and relatively constant in air pressure, temperature and humidity, away from direct sunlight, as well as corrosive, flammable and explosive gases and chemicals. Tilting, vibration, extrusion and other dangerous situations shall be avoided during the placement or handling of the device! Users shall pay attention to the voltage, current and frequency values of the connected power supply, and check whether the grounding terminal is normal. Avoid running the device under overload!

3 - Unpacking and assembly

3.1 Unpacking

Note: The ventilator shall be placed in a well-ventilated environment, away from fire and heat sources, and avoid vibration, direct sunlight and sudden temperature changes. It is recommended to prepare the operating environment according to Section 1.2 - Environment Requirements.

- 1) After unpacking, take out the protective foam and all attachments;
- 2) Please keep all packages and packing materials for future transportation;
- 3) Please check the received materials according to the enclosed packing list. If you find any missing or damaged materials, please contact after-sales personnel of RWD immediately.

3.2 Installation of the bellows components

Precautions before installation:

- Check and confirm that the control unit is in good condition, check the back interface and ensure that the driving gas port and the exhaust port are not blocked;
- Check whether the bellows and bellows cover are in good condition, and make sure that there is no debris inside the bellows.

Note: Before connecting pipes, remove the two rubber plugs shown in Figure 3-1. If no pipes are connected, insert the rubber plugs back to their original positions.


Figure 3-1

Put the folding pouch inside the bellows on the assembly ring of the control unit, gently hold the folding pouch, lift and lower it quickly for several times to straighten the corrugation of the folding pouch. Then cover the bellows cover on the folding pouch, press down gently, rotate the bellows cover clockwise until it is properly stuck in the slot of the base, so as to complete the installation of the bellows and the bellows cover.

Figure 3-2

3.3 Pipeline connection

- Introduction to Y-shaped connecting pipe interface


Figure 3-3

The connection to the R420 ventilator is shown below


Figure 3-4

- 1) As shown in Figure 3-2 and Figure 3-3, connect the B port of the Y-shaped connecting pipe to the 17mm port of the bellows base, connect the A port to the driving gas port of the control unit, and connect the C port to the 30mm port of the bellows;
- 2) Shown as figure 3-4, connect the airway pressure connection tube-Luer-22mm to the pressure transducer port of the control unit, screw it tightly, and then install the 22mm connector (transparent) on the 22mm port of the bellows base;
- 3) Connect the ventilator to the anesthesia system. Remove the breathing bag from the anesthesia machine, connect the 22mm connector (transparent) to the breathing bag interface of the anesthesia machine with a threaded pipe, as shown in Figure 3-5, and close the adjustable pressure limiting valve (APL valve) at the same time;


Figure 3-5

- 4) Connect the ventilator to the gas filter canister. The EXHAUST port of the ventilator can be connected to the gas filter canister through a corrugated pipe.

3.4 Connection of accessories and attachments

Figure 3-6

Name	Description
INSP.HOLD remote control unit	Optional. Used for remote control of INSP.HOLD button
V60 straight joint	Optional. The 15mm port is connected to the middle tube, the 22mm port is connected to the animal lung, and the luer connector is connected to the sampling tube
Ventilator EtCO ₂ watertrap	Optional. Used for dehydration of exhaled air
Ventilator EtCO ₂ sampling tube	Optional. CO ₂ gas collection pipeline

- 1) The INSP.HOLD remote control unit can be used directly after infrared matching, refer to section 4.10.6 for details;

- 2) Connection method of sampling tube, watertrap and V60 straight joint: as shown in Figure 3-6, all interface are marked with letters, among which two ports a are connected to port b and port c respectively, port d is connected to the repeated breathing circuit, and port e is connected to the tracheal intubation.

3.5 Installation method of EtCO₂ watertrap

As shown in Figure 1-1, the EtCO₂ watertrap (optional) is installed on the right side of the ventilator.

The fixing bracket of the watertrap is pre-installed at the factory. When installing, align the watertrap into the installing hole from bottom to top, and then push it horizontally with a little force into the buckle for fixation, as shown in Figure 3-7.

Figure 3-7

4 -Operating instructions

Precautions:

Only the device operating in a specified operating conditions (refer to 1.2-Environment Requirements) can work properly.

4.1 Power-on self-check

The “Self-check” include: turbine function, cooling fan, power supply voltage, battery voltage, flow sensor, pressure sensor, and air tightness.

Preparatory before self-check:

- 1) Adjust the oxygen flow meter on the anesthesia machine to zero;
- 2) Select the circulation loop and close the adjustable pressure limiting valve (APL valve);
- 3) Close animal port of the circulation loop;
- 4) Connect the threaded tube to breathing bag interface of the anesthesia machine;
- 5) Press the quick oxygenation button, and wait until the folding pouch rises to the top of the bellows;
- 6) Press the power switch of the ventilator to automatically enter the self-check interface after initialization (refer to the figure below).

Click [Start] to perform the self-check or click [Home] to skip the self-check and directly enter the main operation interface. It is not recommended to skip the self-check operation.

Note: The interface diagram in this user manual is used for illustration only. The actual parameters are subject to the instrument display.


Figure 4-1

If the self-check is qualified, the interface will prompt “The self-check is normal and will automatically jump to the Home page after 3 seconds!”.

If the self-check result is abnormal, it will be shown as Figure 4-2. The figure shows the detected exception items and detailed information, which is convenient for the users’ troubleshooting. Please refer to section **5-Alarm information and Processing** for processing.

Click [Re-check] to perform the system self-check again, or click [Home] to directly enter the

main interface.


Figure 4-2

4.2 Main interface

Note: since the R420-PRO model product has all the functions of the R420 series ventilator, the subsequent introduction will take the R420-PRO as an example for detailed description.

■ R420


Figure 4-3

■ R420-Plus & R420-Pro




Figure 4-4

List of parameters/icon functions on the main interface

Parameter/icon	Description
	Set alarm items, lower alarm limit and upper alarm limit. Refer to the description in section 4.9 for details
	Click the “Mute” icon to display  , it is in the mute state at this time, and it will be automatically released after 2 minutes of mute. Or click the mute icon again to unmute.
	Click to enter the self-check interface. Refer to the description in section 4.1 for details
	Click to enter the menu interface. The menu item includes 6 sub-items: parameter list, data record, alarm record, Settings, date and time, and manufacturer maintenance.
	Display the battery status
	Displayed only when an external device is connected to the USB interface
	Displayed only when the CO ₂ module requires zero. For details, refer to Section 4.10.4
	Display the current date and time according to the date format and time format in the Settings
	Pressure-time diagram and CO ₂ -time diagram. Refer to section 4.8 for details

-- PEAK cmH2O -- PEEP cmH2O -- VT ml -- MV L/min -- RR bpm -- EtCO2 mmHg -- FiCO2 mmHg	Real-time parameter value display during ventilation
PCV Mode	Display the current ventilation mode. Click to switch VCV/PCV/APNEA modes
15.0 Weight kg	The currently set weight. The system will automatically match other recommended parameters based on the weight input
Tidal volume	Displayed in VCV and APNEA modes. The volume of air per inhalation during mechanical ventilation
RR	Respiratory rate per minute
PEAK	The target pressure in the inspiratory phase or the peak airway pressure in the animal's breathing during mechanical ventilation in PCV
I:E	The ratio of each inspiratory time to expiratory time during mechanical ventilation
P-Trig	In VCV/PCV modes: During mechanical ventilation, in the expiratory phase, when the airway pressure reaches the set trigger pressure, but does not reach the set next inspiratory time, the device will immediately start the assist-control ventilation to enter the inspiratory phase. In APNEA modes: This parameter is used to identify Spontaneous Ventilation or mechanical ventilation
MV	Minute ventilation during the animal's mechanical ventilation
PEEP	During mechanical ventilation, a certain positive pressure is maintained in the airway at the end of the expiration
EtCO ₂	End-tidal CO ₂ concentration of animals
FiCO ₂	Inhalation CO ₂ concentration of animals
APNEA Time	When it is detected that the animal does not breathe or breathes weakly during the [APNEA Time], so as to fail to reach the set value of [P-Trig], APNEA ventilation is initiated immediately

Caution!

In VCV & APNEA mode, there is a linkage relationship among tidal volume, respiratory rate, and inspiration and expiration ratio, and the following conditions must be met: tidal volume (mL) * respiratory rate (bpm) * inspiration and expiration ratio value and (I+E) ≤ 60000, that is, $VT \cdot RR \cdot (I+E) \leq 60000$.

4.3 Ventilation mode

Optional mechanical ventilation modes: VCV mode, PCV mode, APNEA mode.

VCV mode →

15.0 Weight kg	VCV Mode	1500 VT	1:2 I:E	10 PEEP	12 RR	-3 P-Trig	00:00:00 Start
-------------------	-------------	------------	------------	------------	----------	--------------	-------------------

PCV mode →

15.0 Weight kg	PCV Mode	11 PEAK	1:2 I:E	10 PEEP	12 RR	-3 P-Trig	00:00:00 Start
-------------------	-------------	------------	------------	------------	----------	--------------	-------------------

APNEA mode →

15.0 Weight kg	APNEA Mode	1500 VT	1:2 I:E	30 APNEA Time	12 RR	-3 P-Trig	00:00:00 Start
-------------------	---------------	------------	------------	------------------	----------	--------------	-------------------

Figure 4-5

Mode	Introduction
VCV	The mechanical ventilation is triggered by the device or the device ventilation is triggered by the animal's Spontaneous Ventilation, with [tidal volume] as the reference standard and [upper pressure limit] as the overpressure protection.
PCV	The mechanical ventilation is triggered by the device or the device ventilation is triggered by the animal's Spontaneous Ventilation, with [airway peak pressure] as the reference standard.
APNEA	On the basis of the VCV model, the start and stop control of mechanical ventilation is added when the animal appears to breathe on its own.

4.4 Parameter range


After the user selecting the mode on the main interface and setting the animal weight, the system will automatically matches other recommended parameters according to the weight, and the user can also edit the parameters according to actual needs.

The adjustable parameters are shown as follows:

Ventilation parameters	Input range
Weight	0.4~150.0kg
VT (VCV / APNEA mode only)	4~1500 mL
RR	2~60bpm
PEAK (PCV mode only)	5~50cmH ₂ O
I:E	1:1~1:4
P-Trig	1) When PEEP is set to OFF, P-Trig range:-10~ -1cmH ₂ O, OFF 2) When PEEP is set between 3~20 cmH ₂ O, P-Trig range:-10~ -4cmH ₂ O, OFF 3) In APNEA mode, P-Trig range:-10~ -1cmH ₂ O
PEEP (VCV / PCV mode only)	OFF, 3~20 cmH ₂ O In PCV mode, PEAK-PEEP ≥ 5 In VCV mode, PLIMmax-PEEP ≥ 7
APNEA Time (APNEA mode only)	15~30s

4.5 Use of the shuttle knob

Functions of the shuttle knob: rotate to select, click to lock, rotate to edit, and click to confirm.

- 1) Rotate the shuttle knob to select any parameter on the interface, e.g. , rotate clockwise to select right, and rotate counterclockwise to select left.
- 2) To edit the current parameter, first press the shuttle knob to lock the parameter, as shown in the figure , until the setting range appears below the parameter.
- 3) After the parameter is locked, turn the shuttle knob clockwise to increase the parameter value or counterclockwise to decrease the value.
- 4) Press the knob again to confirm the current parameter setting. If there is no confirmation operation within 15s, the current setting will be canceled.
- 5) Rotate the shuttle knob to select the ventilation mode. Firstly rotate the shuttle knob to

select the current mode, for example , click the shuttle knob to pop up the bubble options , then rotate the shuttle knob to select the ventilation mode, click the shuttle knob again to finish the ventilation selection and exit the bubble options.

4.6 Parameter List

Click [Menu]>[Parameter List] on the main interface, as shown in Figure 4-6. This interface will allow users to save and recall common animal breathing parameters. Up to 50 records can be saved.

No. (3/50)	Mode	Weight kg	VT ml	PEAK cmH2O	I:E	RR bpm	PEEP cmH2O	P-Trig cmH2O	APNEA Time s
1	PCV	2.0	--	50	1:1	15	-8	OFF	--
2	VCV	150.0	500	--	1:4	15	-2	OFF	--
3	APNEA	0.4	1500	--	1:3.5	15	--	OFF	15
4	VCV	20.0	343	--	1:2	15	-1	OFF	--
5	PCV	15.0	--	5	1:2	15	-5	OFF	--

Buttons: Add, Edit, Delete, Apply

Figure 4-6

Key name	Key function
	Click to enter the [Add] interface
	Click to enter the [Edit] interface of the selected parameter
	Click to delete the selected parameter
	Click to call the selected parameter and return to the main interface

	Click the button to page up. Long press to quickly page up is allowed
	Display “Current page number/Total pages”
	Click the button to page down. Long press to quickly page down is allowed

4.6.1 Add/edit/delete/call parameters

Click [Add]/[Edit] in the parameter list to enter the interface for adding parameters and editing parameters in Figure 4-7 and Figure 4-8.

In each ventilation mode, the parameter items displayed on the interface are slightly different, and users can switch to select the ventilation mode according to their needs.

The operations for adding and editing parameters are the same. Click / to increase or decrease the value, or select it and use the shuttle knob for locking and editing.

■ Add parameters -VCV


Figure 4-7

■ Edit parameters - PCV


Figure 4-8

For details about how to add or edit parameters in APNEA mode, refer to VCV and PCV modes.

Click [Save] to save the added or edited parameters to the parameter list, or select [Save & Call] to save to the parameter list, and call the adjusted parameters and enter the main interface before start the mechanical ventilation operation.

4.7 Start/end ventilation

After the parameter setting is completed, click  and confirm according to the prompt to start the ventilator ventilation. To end the ventilation, click  to end the ventilation operation.

4.8 Pressure-time diagram & CO₂-time diagram

During ventilation, the interface will display the real-time pressure-time diagram and/or CO₂-time diagram, as shown below.


Figure 4-9

Pressure-time diagram:

The pressure-time diagram will display the real-time pressure value measured by the pressure sensor.

The X-axis range will automatically adjust according to the set respiratory rate.

- When the respiratory rate is set to 26~60bpm, it is 0~10s;
- When the respiratory rate is set to 10~25bpm, it is 0~20s;
- When the respiratory rate is set to 2~9bpm, it is 0~50s;

The Y-axis range will automatically switch the coordinate range with the maximum peak airway pressure shown in the current view..

- When the airway peak pressure is >30 cmH₂O, it is -12~60 cmH₂O;
- When the airway peak pressure is between 16~30 cmH₂O, it is -6~30 cmH₂O;
- When the airway peak pressure is <15 cmH₂O, it is -3~-15 cmH₂O.

CO₂-time diagram:

The CO₂-time diagram will show the exhaled CO₂ value of the animal as measured by the carbon dioxide sensor.

The X-axis range will automatically adjust according to the set respiratory rate.

- When the respiratory rate is set to 26~60bpm, it is 0~10s;
- When the respiratory rate is set to 10~25bpm, it is 0~20s;
- When the respiratory rate is set to 2~9bpm, it is 0~50s;

Y-axis range: 0~80 mmHg.

4.9 Alarm setting

The minute ventilation, upper pressure limit (only VCV mode is valid), and EtCO₂ and FiCO₂ lower alarm limit and upper alarm limit can be set. (The standard ventilator only includes the alarm items of minute ventilation and upper pressure limit.)

Click [Default] to restore the default factory parameters.

Item	Lower	Upper	Default
MV:	16.1 L/min	20.2 L/min	
PLIMmax:		20 cmH2O	
EtCO2:	15 mmHg	20 mmHg	
FiCO2:		20 mmHg	

Figure 4-10

Parameter name	Input range
MV (Not available in the "Spontaneous Ventilation" phase of APNEA mode)	Upper limit: 0.2L/min ~100.0L/min Lower limit: 0.0L/min ~ 99.0L/min Input limits: When the upper limit is ≤ 15.0 L/min, the upper limit value — the lower limit value must be ≥ 0.2 L/min; When the upper limit is > 15.0 L/min, the upper limit value — the lower limit value must be ≥ 1.0 L/min;
Upper pressure limit (Not available in the "Spontaneous Ventilation" phase of APNEA mode, not available in PCV mode)	Upper limit: 2 cmH ₂ O ~ 70 cmH ₂ O
EtCO ₂	Upper limit: OFF, or 2mmHg ~ 99mmHg Input limits: Upper limit value - lower limit value should be ≥ 2 mmHg
	Lower limit: OFF, 0mmHg ~ 97mmHg
FiCO ₂	Upper limit: OFF, 1mmHg ~ 99mmHg

4.10 Settings


Figure 4-11

4.10.1 Brightness

Brightness control: Click the “-” and “+” icons to adjust the system brightness.

4.10.2 Language

Language switch: Click to switch the languages of Chinese/English/Spanish.

4.10.3 Device Zero

Click [Settings] > [Device Zero] to enter the interface shown in the figure below.

Note: Please remove the gas source and remove the breathing system tubes before zeroing. Be sure to unplug the threaded tube connected to the breathing bag interface of the anesthesia machine, refer to Figure 3-5 for instructions.

“Pressure zero” can be performed on the airway pressure sensor or “flow zero” can be performed on the flow sensor. It can be down by clicking the “Clear” button on the right side of the corresponding item. (It is recommended to clear it once a year).

Note: The pressure value must be lower than 2cmH₂O when the airway pressure is cleared. When the flow is cleared, the current flow rate must be >-1L/min and <1L/min.

There will be a prompt in the lower left corner of the corresponding item for successful or unsuccessful clearing. As shown in Figure 4-12, the interface prompts “Airway Pressure cleared successfully!” and “Flow clearing failed!”.

The lower left corner of the interface will display the last clearing time of pressure sensor and flow sensor, which is for user reference only.


Figure 4-12

4.10.4 CO₂ Zero

Click [Settings] > [CO₂ Zero] to enter the interface shown in the figure below.

Note: Before clearing, the device must be turned on for more than 5 minutes, and the EtCO₂ sampling tube must be disconnected.

If it is confirmed to clear CO₂, click “Yes” and wait for the completion of clearing. To exit the current operation, click “No”.


Figure 4-13

4.10.5 CO₂ on/off

Click to **turn on/off** the CO₂ function. It is turned on by default.

Note: If the CO₂ function is turned off, the CO₂ monitoring and related alarms will be turned off.

4.10.6 Remote Pair

Before using the remote control unit, remote pair needs to be completed.

Perform pairing according to the description on the page shown in Figure 4-14. Note this must be completed within 10 seconds:

- 1) Click [Pair].
- 2) Press and hold on the button of the remote controller for 1s.
- 3) If the testspot color changed as button be pressed and released, the match is successful. If the testspot color does not change, the match fails, please repeat step 1~3.
- 4) To cancel the pair, click [Pair Remove], confirm and wait for the pair to be moved.


Figure 4-14

4.11 Usage and remote control of INSP.HOLD button

Note: INSP.HOLD function is not available in APNEA mode.

Hold function is not available in APNEA mode

Scenario: INSP.HOLD function can be used for preoperative induction, CT imaging, thoracic suture, etc.

Button:

- 1) During ventilation, press the INSP.HOLD button on the control unit, the device will maintain the tidal volume (VCV mode)/peak airway pressure (PCV) at the end of the current inspiration, and release the INSP.HOLD;
- 2) Execute the INSP.HOLD function, and display the duration of [Inspiration Hold] and the pressure value of [Airway-Pressure] on the floating window at the end of inspiration;


Figure 4-15

- 3) INSP.HOLD duration: The longest holding time of the INSP.HOLD button is the

INSP.HOLD duration set in [Settings], that is, the time from the start of holding tidal volume/peak airway pressure to the end.

Recommendation: For animals with the weight <3kg, it is recommended that the duration of inspiration hold shall not exceed 5s; for animals with the weight ≥3kg, it is recommended that the duration of inspiration hold shall not exceed 30s.

Precautions:

- 1) When the remote control device is connected, the start/stop function of the INSP.HOLD button can be controlled through remote control;
- 2) During remote control operation, please press the INSP.HOLD button slightly, and do not move left or right or relax the pressing during the pressing, otherwise the connection will be disconnected;
- 3) Do not make the animal's lung expand excessively in the process of the INSP.HOLD taking effect. Excessive inflation of the lung will cause damage to lung organs! It is recommended to preset the appropriate pressure value and tidal volume value before connecting to the animal. It is recommended to use the recommended flow rate and then press and hold the INSP.HOLD button.

4.12 Date & time

Select the current date or time value and enter it through the numeric keyboard on the right to customize the current system date and time.

Click “Date Format” and “Time Format” to display the currently displayed system date and system time according to the selected format.

The screenshot displays the 'Date & Time' configuration screen. At the top, a status bar shows the current date and time as '2021/10/14 11:19 AM'. The main area contains four rows of settings: 'Date' (YYYY-MM-DD) with a value of '2021 / 10 / 14'; 'Time' (H:M) with a value of '11 : 19'; 'Date Format' with three options: 'y/m/d' (selected), 'm/d/y', and 'd/m/y'; and 'Time Format' with two options: '24-Hour System' (selected) and '12-Hour System'. To the right of these settings is a numeric keypad with buttons for digits 1-9, 0, a backspace key (X), and a confirm key (checkmark).

Figure 4-16

5 - Alarm & prompt processing

The alarm information details is as follows:


Figure 5-1

The prompt details is as follows:


Figure 5-2

Alarm signal

In case of an alarm, the alarm signal can be obtained audibly and visually:

- 1) Voice prompt
- 2) Text message

The system distinguish the alarm information level with different background colors:

- 1) High level alarm: red
- 2) Medium level alarm: yellow
- 3) Low level prompt: yellow

Fault code	Alarm level	Alarm information	Solution
001	High level fault	“Animal suffocation”	1) Please check the animal’s vital signs, and restart the device after confirming that the animal is in good condition; 2) Contact after-sales personnel of RWD.
003	High level fault	“Motor stalling!”	1) Self-check again; 2) Restart the device; 3) Please contact the after-sales personnel of RWD.
004	High level fault	“Hall sensor abnormality!”	1) Self-check again; 2) Restart the device; 3) Please contact the after-sales personnel of RWD.
005	High level alarm (R420-PRO only)	“CO ₂ sensor failure!”	1) Restart the device; 2) Please contact the after-sales personnel of RWD.

007	High level alarm	“High airway pressure!”	1) In VCV mode, check the alarm settings: whether [Upper Pressure Limit] is set too low; 2) In VCV mode, check the ventilation settings: whether [Tidal Volume] is set too high; 3) In VCV mode, check whether the animal has Spontaneous Ventilation, whether the animal’s thoracic cavity is compressed, and whether the animal’s lung has changes in compliance; 4) Check the breathing circuit for blockages or kinks; 5) Please restart the mechanical ventilation.
008	High level alarm	“High airway pressure!”	1) Check the breathing circuit for blockages or kinks.
009	High level alarm	“Low airway pressure!”	1) Check the breathing circuit for leaks or disconnections; 2) Check the animal’s trachea cannula for leaks.
101	Medium level alarm	“High ventilation!” minute	1) Check the alarm settings: whether [Upper Minute Ventilation Limit] is set too low; 2) Check the ventilation settings: whether [Airway Peak Pressure], [Tidal Volume], [Respiratory Rate] are set too high; 3) In PCV mode, check the animal’s tracheal intubation for leaks; 4) In PCV mode, check the breathing circuit for leaks.
102	Medium level alarm	“Low ventilation!” minute	1) . Check the alarm settings: whether [Lower Minute Ventilation Limit] is set too high; 2) Check the ventilation settings: whether [Airway Peak Pressure], [Tidal Volume], [Respiratory Rate] are set too low; 3) In PCV mode, check whether the animal’s thoracic cavity is compressed, and whether the animal’s lung has changes in compliance.
103	Medium level alarm (R420-PRO only)	“High EtCO ₂ !”	1) Check the alarm settings: whether [Upper EtCO₂ Limit] is set too low; 2) Check the ventilation settings: whether [Airway Peak Pressure], [Tidal Volume], [Respiratory Rate] are set too low; 3) Check whether the EtCO₂ sampling tube is blocked or kinked; 4) When FiCO₂ is not “0”, check whether the CO₂ absorbent needs to be replaced.
104	Medium level alarm (R420-PRO only)	“Low EtCO ₂ !”	1) Check the alarm settings: whether [Lower EtCO₂ Limit] is set too high; 2) Check the ventilation settings: whether [Airway Peak Pressure], [Tidal Volume], [Respiratory Rate] are set too high;

			3) Check whether the EtCO ₂ sampling tube is disconnected; 4) Check whether the oxygen flow of the anesthesia machine is set too high.
105	Medium level alarm (R420-PRO only)	“High FiCO ₂ !”	1) Check the alarm settings: whether [Upper FiCO ₂ Limit] is set too low; 2) Check the ventilation settings: whether [Airway Peak Pressure], [Tidal Volume], [Respiratory Rate] are set too low; 3) Check whether the CO ₂ absorbent needs to be replaced.
106	Medium level alarm	“Low Battery!”	1) Please connect the power supply.
107	Medium level alarm	“Cooling fan stalling!”	1) Restart the device; 2) Please contact the after-sales personnel of RWD.
108	Medium level alarm	“Power supply error!”	1) Check whether the power supply of ventilator is connected; 2) Restart the device.
109	Medium level alarm	“Abnormal battery voltage!”	1) Check whether the battery power is too low; 2) Check whether the battery connection is abnormal; 3) Restart the device.
110	Medium level alarm (R420-PRO only)	“EtCO ₂ watertrap disconnected”	1) Please connect the EtCO ₂ watertrap.
111	Medium level alarm (R420-PRO only)	“EtCO ₂ sampling tube blocked!”	1) Please check whether the EtCO ₂ sampling tube is abnormal.
112	Medium level alarm (R420-PRO only)	“CO ₂ module overheat!”	1) Please check whether the ambient temperature is too high; 2) Please suspend the use of the CO ₂ function and allow it to cool for 15 minutes before use.
201	Low level prompt (R420-PRO only)	“CO ₂ Module Zero Required!”	Refer to Section 4.10.4 for CO ₂ zero
/	General prompt	“Spontaneous Ventilation”	No processing required
/	General prompt	“Mechanical Ventilation”	No processing required

Caution!

Before the ventilator works normally, it is necessary to have a professional doctor to set the appropriate airway pressure value for the current subject animal. If the tidal volume is too low, the animal will not inhale enough, and if the tidal volume is too high, the animal's chest and lung will be over-inflated, resulting in chest and lung damage.

It is recommended to operate under the guidance of a professional doctor, preset an approximate tidal volume value, observe whether the animal's chest and lung fluctuations, blood oxygen saturation, blood gas and other respiratory parameters are within the ideal range, and make minor adjustments for any deviation!

6 - Troubleshooting

Fault description	Possible cause	Troubleshooting method
When the ventilator is turned on, the power indicator (green light) does not light up	The power is not connected.	1) Plug into a suitable power outlet; 2) Check the output; 3) Check the fuse.
The operation of the ventilator sounds normal, but the indication of the tidal volume is incorrect	There is a problem inside the ventilator.	Contact after-sales personnel of RWD.
The bellows is disengaged from the assembly ring	The lower part of the bellows is damaged or not installed securely.	Use a new bellows or reconnect the bellows.
The ventilator appears to be operating normally, but the bellows cannot be filled up	1) There is a leak of gas in the respiratory system; 2) The gas supply to the anesthesia machine system is insufficient; 3) The bellows is damaged; 4) The bellows is not fully attached to the assembly ring; 5) The safety valve of the anesthesia machine is not fully closed.	1) Check all pipes and piping connections for leaks; 2) Increase the airflow of the anesthesia system; 3) Replace the bellows with a new one; 4) Reconnect the bellows to the assembly ring; 5) Check the safety valve of the anesthesia machine system and close it completely.
Long press the power switch of the membrane panel, and the device cannot be turned on	1) The power cord is not connected; 2) Low battery; 3) Power button of the membrane panel is damaged.	1) Check the connection of the power cord; 2) Make sure that the battery is fully charged. If the battery cannot be used normally, please timely replace it with a new battery; 3) Contact after-sales personnel of RWD.

7 - Cleaning and maintenance

7.1 Product cleaning

Items	Cleaning method	Detailed steps
Surface	Detergent	In the process of cleaning the surface of the ventilator, please wipe the surface of the ventilator housing with a damp cloth ever soaked in a weakly alkaline detergent (water, soapy water with a pH value between 7.0 and 10.5). Caution! The display screen can only be cleaned with a dry, soft and lint-free cloth without liquids.
	Disinfectant	In the process of disinfecting the surface of the ventilator, please wipe the surface of the ventilator housing with a damp cloth ever soaked in a neutral disinfectant (75% alcohol) solution.
	UV light	When the surface of the ventilator is disinfected by UV irradiation, the ventilator shall be placed under a 30W UV lamp for irradiation at a distance of 1m, and the irradiation time shall be shorter than 60min; Caution! The UV radiation is harmful to the human body, please do not stay in the UV room during irradiation.
Pipeline	Clear water	In the process of cleaning the pipeline, please wipe the ventilator pipeline with a damp cloth ever soaked in clear water.

Warning!!!

- In the process of cleaning and sterilizing, please ensure the applicability of the cleaning and sterilization method to each component, as well as the correctness of the cleaning and sterilization method.
- All liquids shall be kept away from electronic components. Do not allow liquids to penetrate into the device housing.
- Make sure the relevant parts after cleaning or disinfection are completely dry before reconnecting the AC power.
- Do not use rough cloths or scented cleaners.
- Do not use abrasive cleaners (such as steel wool, silver polish or detergent). The PH of the cleaning solution must be 7.0 ~ 10.5.
- Do not sterilize the bellows cover with high-temperature steam, which will result in deformation and failure of normal use.
- In the process of cleaning and disinfecting the circuit, do not wipe the wire terminals directly, otherwise the wire terminals may be oxidized or corroded.
- The folding pouch can only be cleaned with clear water.

7.2 Battery performance maintenance

During use, the battery should be regularly optimized to maintain its lifespan. It is recommended to perform battery maintenance every 3 months of storage.

Please follow the steps below:

- 1) Take out the device and keep it stable for more than 2H at room temperature of 5°C~30°C.
- 2) Switch on the external power supply and charge the battery for 5~7H.
- 3) Reposition the device.
- 4) Battery capacity: 5,000 mAh.

Check the battery performance

The battery performance may degrade over time. Under normal circumstances, the battery with a full charge can be used for 2 hours. If the battery service time is significantly less than 2 hours, users shall consider replacing the battery or contact after-sales personnel of RWD.

Note: The service life of battery depends on the frequency and time of use. The lithium battery properly maintained and stored can be used for about 2 years. Lifespan of the battery not being used improperly may be shorter. We recommend to replace the lithium battery every 2 years.

Battery recycling

If the battery is obviously damaged, or the battery capacity is exhausted, it should be replaced and recycled properly. Appropriate regulations should be followed for disposing of used batteries.

Warning!!!

Do not remove the battery, put it into a fire, or short-circuit it. Battery burning, explosion and leakage may cause personal injury.

7.3 Device zero and CO₂ zero

Refer to the descriptions in the Sections 4.10.3 and 4.10.4.

7.4 Battery replacement

- 1) Unscrew the battery compartment on the back of the machine with a cross screwdriver, as shown in Figure 7-1;


Figure 7-1

Figure 7-2

- 2) Remove the battery cover and put it aside with the screws to observe the wiring harness of the battery; pull out the wire clip shown in Figure 7-2 out of the cabinet;
- 3) Slowly pull out the lithium battery pack (**be careful to prevent the wire harness from getting stuck in the process of pulling, resulting in the rupture of the insulation layer, thereby damaging the battery pack**);
- 4) After pulling out the whole lithium battery pack, loosen the connection of the wire clip shown in the figure, and replace with a new lithium battery pack;
- 5) Finally, reverse the above steps, put the whole new lithium battery package into the battery compartment of the cabinet, and cover the battery cover.

7.5 Fuse replacement

Before replacing the fuse, please switch off and unplug the power plug! Only the same type of fuse can be used, otherwise the instrument may be seriously damaged.

Replacement method: pry open the fuse box on the back of the instrument with a tweezer along the buckle position indicated by the arrow, take out the old fuse, and replace it with a new one.


Figure 7-3

Figure 7-4

7.6 Replacement of turbine air filter

Note: It is recommended to replace the turbine filter once a year, and check and clean the first effect cotton every three months.

- 1) First, loosen and remove the screw counterclockwise at the position of the air filter on the back of the ventilator (Figure 7-1). If the screw is too tight, remove it with a screwdriver, and open the cover, as shown in Figure 7-2;

Figure 7-5

Figure 7-6

- 2) Open the buckle and take out the entire turbine filter.

Figure 7-7


Figure 7-8

7.7 Replacement of EtCO₂ watertrap and EtCO₂ extension tube

Figure 7-9

Figure 7-10

8 - Warranty

RWD provides after-sales support such as device maintenance and parts replacement in case of the device failure due to defects in materials or process during the warranty period.

Any device damage caused by incorrect use or out-of-range use is out of the scope of the warranty, and in case of any need for repairs or parts replacement, the costs incurred should be borne by the user.

If the device returned for repair is found to have been disassembled without authorization of RWD upon arrival, RWD will not provide after-sales services such as quality warranty, free maintenance and parts replacement.

The warranty statement (including its limitations) is issued exclusively by RWD and covers all other warranty conditions.

Appendix I

Intelligent parameters - alarm items and default values

Alarm item	Default value
Lower EtCO ₂ limit mmHg	25
Upper EtCO ₂ limit mmHg	50
Upper FiCO ₂ limit mmHg	4
Upper pressure limit cmH ₂ O	Vary with the input weight, refer to the table below
Upper minute ventilation limit L/min	Vary with the input weight, refer to the table below
Lower minute ventilation limit L/min	Vary with the input weight, refer to the table below

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
0.4	13	0.4	0.0
0.5	13	0.5	0.1
0.6	13	0.6	0.1
0.7	13	0.7	0.1
0.8	13	0.8	0.1
0.9	13	0.9	0.1
1.0	13	1.0	0.1
1.1	13	1.1	0.1
1.2	13	1.2	0.1
1.3	13	1.3	0.1
1.4	13	1.4	0.1
1.5	13	1.5	0.2
1.6	13	1.6	0.2
1.7	13	1.8	0.2
1.8	13	1.9	0.2
1.9	13	2.0	0.2
2.0	13	2.1	0.2
2.1	13	2.2	0.2
2.2	13	2.3	0.2
2.3	13	2.4	0.2
2.4	13	2.5	0.2
2.5	13	2.6	0.3
2.6	13	2.7	0.3
2.7	13	2.8	0.3
2.8	13	2.9	0.3
2.9	13	3.0	0.3
3.0	13	3.1	0.3

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
3.1	15	3.2	0.3
3.2	15	3.3	0.3
3.3	15	3.4	0.3
3.4	15	3.5	0.4
3.5	15	3.6	0.4
3.6	15	3.7	0.4
3.7	15	3.8	0.4
3.8	15	3.9	0.4
3.9	15	4.0	0.4
4.0	15	4.1	0.4
4.1	15	4.2	0.4
4.2	15	4.3	0.4
4.3	15	4.4	0.4
4.4	15	4.5	0.5
4.5	15	4.6	0.5
4.6	15	4.7	0.5
4.7	15	4.8	0.5
4.8	15	4.9	0.5
4.9	15	5.0	0.5
5.0	15	5.1	0.5
5.1	15	3.1	0.5
5.2	15	3.1	0.5
5.3	15	3.2	0.5
5.4	15	3.2	0.5
5.5	15	3.3	0.5
5.6	15	3.2	0.5
5.7	15	3.2	0.5
5.8	15	3.3	0.5
5.9	15	3.3	0.5
6.0	15	3.4	0.5
6.1	15	3.2	0.5
6.2	15	3.2	0.5
6.3	15	3.3	0.5
6.4	15	3.3	0.5
6.5	15	3.4	0.5
6.6	15	3.2	0.5
6.7	15	3.3	0.5
6.8	15	3.3	0.5
6.9	15	3.3	0.5
7.0	15	3.4	0.5

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
7.1	15	3.2	0.5
7.2	15	3.2	0.5
7.3	15	3.3	0.5
7.4	15	3.3	0.5
7.5	15	3.4	0.5
7.6	15	3.1	0.5
7.7	15	3.2	0.5
7.8	15	3.2	0.5
7.9	15	3.3	0.5
8.0	15	3.3	0.5
8.1	15	3.3	0.5
8.2	15	3.4	0.5
8.3	15	3.4	0.5
8.4	15	3.5	0.6
8.5	15	3.5	0.6
8.6	15	3.5	0.6
8.7	15	3.6	0.6
8.8	15	3.6	0.6
8.9	15	3.7	0.6
9.0	15	3.7	0.6
9.1	15	3.7	0.6
9.2	15	3.8	0.6
9.3	15	3.8	0.6
9.4	15	3.8	0.6
9.5	15	3.9	0.6
9.6	15	3.9	0.6
9.7	15	4.0	0.6
9.8	15	4.0	0.6
9.9	15	4.0	0.6
10.0	15	4.1	0.7
10.5	16	4.1	0.7
11.0	16	4.3	0.7
11.5	16	4.5	0.7
12.0	16	4.7	0.7
12.5	16	4.8	0.8
13.0	16	5.0	0.8
13.5	16	5.2	0.8
14.0	16	5.4	0.9
14.5	16	5.6	0.9
15.0	16	5.7	0.9

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
15.5	16	5.7	0.9
16.0	16	5.9	0.9
16.5	16	6.1	1.0
17.0	16	6.2	1.0
17.5	16	6.4	1.0
18.0	16	6.6	1.1
18.5	16	6.7	1.1
19.0	16	6.9	1.1
19.5	16	7.1	1.1
20.0	16	7.2	1.2
20.5	20	7.3	1.2
21.0	20	7.4	1.2
21.5	20	7.6	1.2
22.0	20	7.7	1.2
22.5	20	7.9	1.3
23.0	20	8.0	1.3
23.5	20	8.2	1.3
24.0	20	8.3	1.3
24.5	20	8.5	1.4
25.0	20	8.6	1.4
25.5	20	8.8	1.4
26.0	20	8.9	1.4
26.5	20	9.1	1.4
27.0	20	9.2	1.5
27.5	20	9.4	1.5
28.0	20	9.5	1.5
28.5	20	9.7	1.5
29.0	20	9.8	1.6
29.5	20	10.0	1.6
30.0	20	10.1	1.6
30.5	20	7.7	1.2
31.0	20	7.8	1.2
31.5	20	7.9	1.3
32.0	20	8.0	1.3
32.5	20	8.1	1.3
33.0	20	8.3	1.3
33.5	20	8.4	1.3
34.0	20	8.5	1.4
34.5	20	8.6	1.4
35.0	20	8.7	1.4

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
35.5	23	8.8	1.4
36.0	23	8.9	1.4
36.5	23	9.0	1.4
37.0	23	9.2	1.5
37.5	23	9.3	1.5
38.0	23	9.4	1.5
38.5	23	9.5	1.5
39.0	23	9.6	1.5
39.5	23	9.7	1.6
40.0	23	9.8	1.6
40.5	23	9.9	1.6
41.0	23	10.1	1.6
41.5	23	10.2	1.6
42.0	23	10.3	1.6
42.5	23	10.4	1.7
43.0	23	10.5	1.7
43.5	23	10.6	1.7
44.0	23	10.7	1.7
44.5	23	10.8	1.7
45.0	23	11.0	1.8
45.5	23	11.1	1.8
46.0	23	11.2	1.8
46.5	23	11.3	1.8
47.0	23	11.4	1.8
47.5	23	11.5	1.8
48.0	23	11.6	1.9
48.5	23	11.7	1.9
49.0	23	11.9	1.9
49.5	23	12.0	1.9
50.0	23	12.1	1.9
50.5	23	12.2	1.2
51.0	23	12.3	1.2
51.5	23	12.4	1.2
52.0	23	12.5	1.2
52.5	23	12.6	1.3
53.0	23	12.7	1.3
53.5	25	12.8	1.3
54.0	25	12.9	1.3
54.5	25	13.0	1.3
55.0	25	13.1	1.3

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
55.5	25	13.2	1.3
56.0	25	13.3	1.3
56.5	25	13.4	1.3
57.0	25	13.5	1.4
57.5	25	13.6	1.4
58.0	25	13.7	1.4
58.5	25	13.8	1.4
59.0	25	13.9	1.4
59.5	25	14.0	1.4
60.0	25	14.1	1.4
60.5	28	14.2	1.4
61.0	28	14.3	1.4
61.5	28	14.4	1.4
62.0	28	14.5	1.5
62.5	28	14.6	1.5
63.0	28	14.7	1.5
63.5	28	14.8	1.5
64.0	28	14.9	1.5
64.5	28	15.0	1.5
65.0	28	15.1	1.5
65.5	28	15.2	1.5
66.0	28	15.3	1.5
66.5	28	15.4	1.5
67.0	28	15.5	1.6
67.5	28	15.6	1.6
68.0	28	15.7	1.6
68.5	28	15.8	1.6
69.0	28	15.9	1.6
69.5	28	16.0	1.6
70.0	28	16.1	1.6
70.5	28	16.2	1.6
71.0	28	16.3	1.6
71.5	28	16.4	1.6
72.0	28	16.5	1.7
72.5	28	16.6	1.7
73.0	28	16.7	1.7
73.5	28	16.8	1.7
74.0	28	16.9	1.7
74.5	28	17.0	1.7
75.0	28	17.1	1.7

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
75.5	28	17.2	1.7
76.0	28	17.3	1.7
76.5	28	17.4	1.7
77.0	28	17.6	1.8
77.5	28	17.7	1.8
78.0	28	17.8	1.8
78.5	28	17.9	1.8
79.0	28	18.0	1.8
79.5	28	18.1	1.8
80.0	28	18.2	1.8
80.5	28	18.3	1.8
81.0	28	18.4	1.8
81.5	28	18.5	1.8
82.0	28	18.6	1.9
82.5	28	18.7	1.9
83.0	28	18.8	1.9
83.5	28	18.9	1.9
84.0	28	19.0	1.9
84.5	28	19.1	1.9
85.0	28	19.2	1.9
85.5	28	19.3	1.9
86.0	28	19.4	1.9
86.5	28	19.5	1.9
87.0	28	19.6	2.0
87.5	28	19.7	2.0
88.0	28	19.8	2.0
88.5	28	19.9	2.0
89.0	28	20.0	2.0
89.5	28	20.1	2.0
90.0	28	20.2	2.0
90.5	28	20.3	2.0
91.0	28	20.4	2.0
91.5	28	20.5	2.0
92.0	28	20.6	2.1
92.5	28	20.7	2.1
93.0	28	20.8	2.1
93.5	28	20.9	2.1
94.0	28	21.0	2.1
94.5	28	21.1	2.1
95.0	28	21.2	2.1

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
95.5	28	21.3	2.1
96.0	28	21.4	2.1
96.5	28	21.5	2.1
97.0	28	21.6	2.2
97.5	28	21.7	2.2
98.0	28	21.8	2.2
98.5	28	21.9	2.2
99.0	28	22.0	2.2
99.5	28	22.1	2.2
100.0	28	22.2	2.2
100.5	28	22.3	2.2
101.0	28	22.4	2.2
101.5	28	22.5	2.3
102.0	28	22.6	2.3
102.5	28	22.7	2.3
103.0	28	22.8	2.3
103.5	28	22.9	2.3
104.0	28	23.0	2.3
104.5	28	23.1	2.3
105.0	28	23.2	2.3
105.5	28	23.3	2.3
106.0	28	23.4	2.3
106.5	28	23.5	2.4
107.0	28	23.6	2.4
107.5	28	23.7	2.4
108.0	28	23.8	2.4
108.5	28	23.9	2.4
109.0	28	24.0	2.4
109.5	28	24.1	2.4
110.0	28	24.2	2.4
110.5	28	24.3	2.4
111.0	28	24.4	2.4
111.5	28	24.5	2.5
112.0	28	24.6	2.5
112.5	28	24.7	2.5
113.0	28	24.8	2.5
113.5	28	24.9	2.5
114.0	28	25.0	2.5
114.5	28	25.1	2.5
115.0	28	25.2	2.5

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
115.5	28	25.3	2.5
116.0	28	25.4	2.5
116.5	28	25.5	2.6
117.0	28	25.7	2.6
117.5	28	25.8	2.6
118.0	28	25.9	2.6
118.5	28	26.0	2.6
119.0	28	26.1	2.6
119.5	28	26.2	2.6
120.0	28	26.3	2.6
120.5	28	26.4	2.6
121.0	28	26.5	2.6
121.5	28	26.6	2.7
122.0	28	26.7	2.7
122.5	28	26.8	2.7
123.0	28	26.9	2.7
123.5	28	27.0	2.7
124.0	28	27.1	2.7
124.5	28	27.2	2.7
125.0	28	27.3	2.7
125.5	28	27.4	2.7
126.0	28	27.5	2.7
126.5	28	27.6	2.8
127.0	28	27.7	2.8
127.5	28	27.8	2.8
128.0	28	27.9	2.8
128.5	28	28.0	2.8
129.0	28	28.1	2.8
129.5	28	28.2	2.8
130.0	28	28.3	2.8
130.5	28	28.4	2.8
131.0	28	28.5	2.8
131.5	28	28.6	2.9
132.0	28	28.7	2.9
132.5	28	28.8	2.9
133.0	28	28.9	2.9
133.5	28	29.0	2.9
134.0	28	29.1	2.9
134.5	28	29.2	2.9
135.0	28	29.3	2.9

Weight kg	Upper pressure limit cmH ₂ O	Upper minute ventilation limit L/min	Lower minute ventilation limit L/min
135.5	28	29.4	2.9
136.0	28	29.5	2.9
136.5	28	29.6	3.0
137.0	28	29.7	3.0
137.5	28	29.8	3.0
138.0	28	29.9	3.0
138.5	28	30.0	3.0
139.0	28	30.1	3.0
139.5	28	30.2	3.0
140.0	28	30.3	3.0
140.5	28	30.4	3.0
141.0	28	30.5	3.1
141.5	28	30.6	3.1
142.0	28	30.7	3.1
142.5	28	30.8	3.1
143.0	28	30.9	3.1
143.5	28	31.0	3.1
144.0	28	31.1	3.1
144.5	28	31.2	3.1
145.0	28	31.3	3.1
145.5	28	31.4	3.1
146.0	28	31.5	3.2
146.5	28	31.6	3.2
147.0	28	31.7	3.2
147.5	28	31.8	3.2
148.0	28	31.9	3.2
148.5	28	32.0	3.2
149.0	28	32.1	3.2
149.5	28	32.2	3.2
150.0	28	32.3	3.2

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The screenshot displays the 'Fecha y hora' (Date and Time) configuration interface. At the top, there is a back arrow, the title 'Fecha y hora', and a status bar showing the current date and time '2021/10/14 11:19 AM' along with a battery icon. The main area contains four rows of settings:

- Fecha** (YYYY-MM-DD): A text field showing '2021 / 10 / 14'.
- Hora** (H:M): A text field showing '11 : 19'.
- Formato de fecha**: Three buttons labeled 'y/m/d', 'm/d/y', and 'd/m/y'. The 'd/m/y' button is currently selected.
- Formato de hora**: Two buttons labeled '24-Hour System' and '12-Hour System'. The '24-Hour System' button is currently selected.

To the right of these settings is a numeric keypad with buttons for digits 1-9, 0, a backspace key (X), and a confirmation key (checkmark).

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Share Info GmbH