

AGSS Gas Supply



Waste gas scavenging

Three types of waste gas scavenging system, fit for different scenarios.

➤ **Active AGSS**
Suitable for hospitals
with active waste gas
suction system




➤ **Passive AGSS**
Exhaust waste gas
from OR through the
tube



➤ **absorber**
Trolley bracket with tray



Waste gas scavenging

Pic	Part number	Description	Apply to
	115-002342-00	Passive AGSS accessory kit	All

Length: 4.7m



including:

-Tube, 32mm, 1 pcs

-Adaptor, 1 pcs



Waste gas scavenging

Pic	Part number	Description	Apply to
	115-011860-00	AGSS kit, low flow	WATO EX-20Vet
	115-011859-00	AGSS kit, high flow	
	115-030332-00	AGSS kit, low flow	WATO EX-35Vet WATO EX-65Pro Vet
	115-030333-00	AGSS kit, high flow	
	115-076014-00	AGSS kit, low flow	Veta 5/3
	115-076013-00	AGSS kit, high flow	

Waste gas scavenging

Active AGSS

AGSS transfer and receiving system	
Size	430 mm×132 mm×114 mm (height x width x depth)
Type of disposal system	High-flow disposal system
Pump rate	75 to 105 L/min

AGSS transfer and receiving system	
Size	430 mm×132 mm×114 mm (height x width x depth)
Type of disposal system	Low-flow disposal system
Pump rate	25 to 50 L/min

AGSS transfer hose assembly	
Length	0.3 m
AGSS active drain tube assembly	
Length	4 m

