

D-Pro&D-U2000 veterinary dental unit was designed to help veterinarians perform their daily work easier and more efficient. Veterinarians can complete grinding, polishing, cutting and cleaning very quickly without fatigue.

With RWD dental unit, you will...

- Get high return on equipment investment.
- Shorten the treatment, therefore reduce the risk of anesthesia.
- Reduce the damage of alveolar bone and oral mucosa caused by violent tooth extraction.
- Reduce the risk of retained roots.
- Make more flexible dental treatment solutions.

D-Pro Veterinary Dental Unit >>>



360° handpiece, light in weight, reduces hand fatigue effectively.



Low speed motor with external water for effective cooling. It can also be compatible with polishing head and burs.



High-speed LED handpiece provides bright operative vision.



Automatically drain pollution or water from the air tank, greatly reduces the trouble of daily maintenance, reduces the risk of rusting, and extends the lifespan.



20L large capacity air tank ensures pressure balance and long-term supply, and reduces the frequency of air charging.



High quality oil-free air compressor to avoid the trouble of daily maintenance.



Magnetic holder to prevent fine dental tools from accidental loss.




The whole machine weighs only 46kg, flexible movement with super smooth casters.




Ultrasonic Scaler Coolant Outlet




High Speed Handpiece



Magnetic Holder



Low Speed Polisher



3 Way Air/Water Syringe



Veterinary Ultrasonic Scaler (Optional)