


D-U2000 Flexible Veterinary Dental Unit >>>



Powerful 360° swivel high-speed handpiece with self-generating LED illumination helps veterinarians cut animal teeth easily.



Automatically drain pollution or water from the air tank, greatly reduces the trouble of daily maintenance, reduces the risk of rusting, and extends the lifespan.



Low speed motor with coolant delivery tube for effective cooling. It can also compatible with polishing head and burs.



Magnetic platform to absorb burs or other fine tools, avoid accidental sliding or loss.



28000Hz scaler with automatic frequency tracking to ensure power, and white LED light for clear vision without fatigue.



Whole machine weighs only 32kg. Super smooth casters will ensure flexible movement.




High-speed handpiece with white light



Low-speed handpiece with coolant delivery tube



Operation console



3 way air/ water syringe



Magnetic platform



Powerful scaler